BMG

Immuno Enzymometric Assay

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Attention

For North and South American Customers: Please refer to the AIA-AAM Docs on CD for the appropriate information.

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ST AIA-PACK BMG

For Quantitative Measurement of β_2 microglobulin (BMG) in Serum, Heparinized Plasma or Urine

NAME AND INTENDED USE

ST AIA-PACK BMG is designed for IN VITRO DIAGNOSTIC USE ONLY for the quantitative measurement of β_2 microglobulin (BMG) in human serum, heparinized plasma or urine on TOSOH AIA System Analyzers.

SUMMARY AND EXPLANATION OF TEST

Human β_2 microglobulin (BMG) is a single polypeptide of 100 amino acids, identified in 1968 by Berggard et al (1). The molecular weight of BMG is 11,800 daltons and represents the light-chain moiety of the major histocompatibility complex class antigens (2-4). It is worth noting that BMG shows approximately 35% sequence homology with immunoglobulin domains. BMG is shed from the cell surface of all nucleated cells (2,3) during growth and differentiation, and cleared in the kidneys by filtration through the glomeruli. Since it is reabsorbed by the proximal tubular cells (5) and degraded, only trace amounts are detected in urine from healthy individuals, whereas marked elevation is seen in patients with proximal tubular dysfunctions due to exposure to aminoglycosides (6), anti-cancer or inflammatory drugs (7,8) and heavy metals (9). Determination of serum BMG is of diagnostic value in a variety of disorders. It can be used to assess glomerular function in advanced renal disease (10,11) and after renal transplantation (12,13). It can provide a measure of tumor burden and prognosis in patients with multiple myeloma (14,15), B cell lymphoma (16) and chronic lymphocytic leukemia (17). It is useful to monitor disease activity in chronic inflammatory disorders, such as rheumatoid arthritis (18), systemic lupus erythematosus (19), Sjogren's syndrome (20), and Crohn's disease (21). More recently, BMG has been reported to be of value in monitoring patients with HIV infection (22-25).

PRINCIPLE OF THE ASSAY

The ST AIA-PACK BMG is a two-site immunoenzymometric assay which is performed entirely in the AIA-PACK test cups. β_2 microglobulin present in the test sample is bound with monoclonal antibody immobilized on a magnetic solid phase and enzyme-labeled polyclonal antibody in the AIA-PACK test cups. The magnetic beads are washed to remove unbound enzyme-labeled monoclonal antibody and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labeled monoclonal antibody that binds to the beads is directly proportional to the β_2 microglobulin concentration in the test sample. A standard curve is constructed, and unknown sample concentrations are calculated using this curve.

MATERIAL PROVIDED (ST AIA-PACK BMG, Cat. No. 0025259)

Plastic test cups containing lyophilized twelve magnetic beads coated with anti- β_2 microglobulin mouse monoclonal antibody and 135 μ L of anti- β_2 microglobulin rabbit polyclonal antibody conjugated to bovine alkaline phosphatase with sodium azide as a preservative.

MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are not provided but are required to perform β_2 microglobulin analysis using the ST AIA-PACK BMG (Cat. No. 0025259) on the TOSOH AIA System Analyzers. They are available separately from TOSOH.

Materials	Cat. No.
AIA Nex-IA or AIA-21	0018539
AIA Nex-IA or AIA-21 LA	0018540
AIA-1800 ST	0019836
AIA-1800 LA	0019837
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-600 II	0019014
AIA-600 II BCR	0019328
AIA-360	0019945
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK SUBSTRATE REAGENT II	
AIA-PACK SUBSTRATE RECONSTITUENT II	
AIA-PACK BMG CALIBRATOR SET	0020359
CALIBRATOR (1)	0 mg/L
CALIBRATOR (2)	0.0125 mg/L (approx.)
CALIBRATOR (3)	0.05 mg/L (approx.)
CALIBRATOR (4)	0.15 mg/L (approx.)
CALIBRATOR (5)	0.30 mg/L (approx.)
CALIBRATOR (6)	0.42 mg/L (approx.)
AIA-PACK BMG SAMPLE DILUTING CONCENTRATE	0020559
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
SAMPLE CUPS	0018581
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971

Additional Requirements for AIA Nex-IA / AIA-21 only:

PIPETTE TIPS	0018552
PRELOADED PIPETTE TIPS	0018583
Additional Requirements for AIA-600 II, AIA-1800 and AIA-2000:	
PIPETTE TIPS	0019215
TIP RACK	0019216
PRELOADED PIPETTE TIPS	0022103

Only materials obtained from TOSOH should be used. Materials obtained elsewhere should not be substituted since assay performance is characterized based strictly on TOSOH materials.

WARNINGS AND PRECAUTIONS

- The ST AIA-PACK BMG is intended for in vitro diagnostic use only.
- Test cups from different lots or different assays should not be mixed within a tray.
- The ST AIA-PACK BMG contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human serum is not used in the preparation of this product; however, since human specimens will be used for samples and other quality control products in the lab may be derived from human serum, please use standard laboratory safety procedures in handling all specimens and controls.
- Do not use beyond the expiration date.
- The ST AIA-PACK BMG has been designed so that the high dose "hook effect" is not a problem for the vast majority of samples. On the recommended 51-fold dilution, serum or heparinized plasma samples with β_2 microglobulin concentrations between 20 and 102 mg/L will read > 20 mg/L. On the recommended 5-fold dilution, urine samples with β_2 microglobulin concentration between 2.0 and 10.0 mg/L will read > 2.0 mg/L. Using the recommended initial dilutions, the "hook effect" phenomenon may occur at β_2 microglobulin concentrations > 102 mg/L in serum or heparinized plasma, and > 10.0 mg/L in urine.
- TOSOH AIA-1200 series or AIA-600 Immunoassay Analyzers can NOT be used to perform the ST AIA-PACK BMG assay.

STORAGE AND STABILITY

All unopened materials are stable until the expiration date on the label when stored at the specified temperature.

Materials	Cat. No.
2-8°C:	
ST AIA-PACK BMG	0025259
AIA-PACK BMG CALIBRATOR SET	0020359
AIA-PACK BMG SAMPLE DILUTING CONCENTRATE	0020559
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
1-30°C:	
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971

ST AIA-PACK BMG test cups may be stored at 18-25°C for up to 1 day. Calibrators and Sample Diluting Concentrate should be used within 1 day and 7 days of opening or reconstituting, respectively, provided the vials or bottles are kept tightly sealed and refrigerated at 2-8°C. They should be used after equilibrating to 18-25°C for about 30 minutes. The BMG Sample Diluting Solution prepared following the instructions in the I. Reagent Preparation, PROCEDURE section below can be used for up to 90 days provided that 1) it is used with only ST AIA-PACK BMG, 2) it is used for manual dilutions ONLY, and 3) the bottles are kept tightly sealed and refrigerated immediately after use. Reconstituted substrate solution is stable for 3 days at 18-25°C or 7 days at 2-8°C. Working diluent and wash solutions are stable for 30 days at 18-25°C. Reagents should not be used if they appear cloudy or discolored.

SPECIMEN COLLECTION AND HANDLING

- Serum or heparinized plasma or urine is required for the assay. EDTA and citrated plasma SHOULD NOT BE USED.
- If using serum, a venous blood sample is collected aseptically without additives. Store at 18-25°C until a clot has formed (usually 15-45 minutes), then centrifuge to obtain the serum specimen for assay.
- If using heparinized plasma, a venous blood sample is collected aseptically with designated additive. Centrifuge and separate plasma from the packed cells as soon as possible.
- Urine should be collected in a container to which no preservatives have been added. After collection the urine pH should be adjusted to read 7-8 if the testing will not be done immediately.
- Specimen types should not be used interchangeably during serial monitoring of an individual patient. Measured concentrations may vary slightly between sample types in certain patients.
- Samples may be stored at 2-8°C for up to 24 hours prior to analysis. If the analysis cannot be done within 24 hours, the sample should be stored frozen at -20°C or below for up to 60 days.
- Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged prior to testing. Prior to assay, bring frozen samples to 18-25°C slowly and mix gently.
- All patient specimens require dilution prior to analysis. Serum or heparinized plasma specimens are diluted by 51 fold and urine specimens by 5 fold with the BMG Sample Diluting Solution prepared following the instructions in the I. Reagent Preparation, PROCEDURE section below. The dilution factors are identified in the test file as Sp.1 for serum or heparinized plasma (51) and Sp.2 for urine (5). The AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800 and AIA-2000 will automatically perform these dilutions and calculate the results. For further information regarding instrument operation, consult the specific TOSOH AIA System Operator's Manual. AIA-PACK BMG CALIBRATOR materials do not require dilution prior to assay.
- The diluted sample required for analysis is 15 μ L of serum or heparinized plasma or urine.

PROCEDURE

For the AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800, AIA-2000 and AIA-360, please refer to their Operator's Manual for detailed instructions.

I. Reagent Preparation

A) Substrate Solution

Bring all reagents to 18-25°C before preparing the working reagent. Add the entire contents of the AIA-PACK SUBSTRATE RECONSTITUENT II (100 mL) to the lyophilized AIA-PACK SUBSTRATE REAGENT II and mix thoroughly to dissolve the solid material.

B) Wash Solution

Add the entire contents of the AIA-PACK WASH CONCENTRATE (100 mL) to approximately 2.0 L of CAP Class I or NCCLS Type I Reagent Grade water, mix well, and adjust the final volume to 2.5 L.

C) Diluent

Add the entire contents of the AIA-PACK DILUENT CONCENTRATE (100 mL) to approximately 4.0 L of CAP Class I or NCCLS Type I Reagent Grade water, mix well, and adjust the final volume to 5.0 L.

D) BMG Sample Diluting Solution

Add 90 mL of CAP Class I or NCCLS Type I Reagent Grade water to the AIA-PACK BMG SAMPLE DILUTING CONCENTRATE (10 mL) and mix well.

II. Calibration Procedure

A) Calibration Curve

The calibrators for use with the ST AIA-PACK BMG are prepared gravimetrically and are compared to internal reference standards.


The calibration curve for ST AIA-PACK BMG is stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and TOSOH AIA System maintenance according to the manufacturer's instructions.

Recalibration may be necessary more frequently if controls are out of the established range for this assay or when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustment or change). For further information regarding instrument operation, consult the TOSOH AIA System Operator's Manual.

A sample calibration curve from AIA-1800 follows and shows the algorithm used for calculating results. The concentration range of the calibration curve is displayed in 1/51 of the assay range in serum or heparinized plasma specimens and in 1/5 of assay range in urine specimens. The concentrations of patient specimens are calculated by multiplying the concentrations obtained on the calibration curve with the dilution factor. AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800 and AIA-2000 will automatically calculate the concentrations of patient specimens using the dilution factor and report the results.

[#BMG-034] Calibrator Lot: 2000405406070000

$$Y = Ax^2 + Bx^2 + Cx + D$$

B) Calibration Procedure

1. Refer to the appropriate TOSOH AIA System Operator's Manual for the procedural instructions.
2. Verify that both the calibrator lot and concentration numbers have been correctly entered into the software.
3. CALIBRATOR (1) for ST AIA-PACK BMG is provided ready for use. TOSOH recommends that all calibrators be run in triplicate.
4. CALIBRATOR (2)-(6) for ST AIA-PACK BMG are lyophilized. All levels should be reconstituted with 1.0 mL of CAP Class I or NCCLS Type I Reagent Grade water. TOSOH recommends that all calibrators be run in triplicate.

C) Calibration Acceptability Criteria

1. The mean rate for the CALIBRATOR (1) should be < 3.0 nmol/(L·s).
2. Since there is a direct relationship between concentration and rate, the rate should increase as the concentration increases.
3. The replicate values should be within a 10% range.

D) Calibration Review and Acceptance

1. Review the calibration curve carefully, using the criteria listed above.
2. Edit the calibration if necessary, then accept the calibration.

For further information regarding calibration, consult the TOSOH AIA System Operator's Manual.

III. Quality Control Procedure

A) Commercially Available Controls

Commercially available controls should be run at least once per day. It is recommended that at least two levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay designates the following:

Control Material: _____

Frequency: _____

Lot number of control material, acceptable limits, and corrective action to be taken if controls do not meet laboratory criteria will be found in a separate quality control document maintained by the laboratory.

B) Quality Control Procedure

1. Assay quality control specimens as instructed in the specific Operator's Manual for your analyzer. In addition, refer to the TOSOH AIA System Operator's Manual for detailed instructions on defining and editing the files.
2. Quality control material to be run with this assay is defined by individual laboratory policy.

IV. Specimen Processing

A) Preparation

All patient specimens require dilution prior to analysis. Serum or heparinized plasma specimens are diluted by 51 fold and urine specimens by 5 fold with the BMG Sample Diluting Solution. Either diluted or undiluted patient specimens are placed in appropriate positions on the TOSOH AIA System Analyzers. For further information regarding sample preparation and instrument operation, refer to the SPECIMEN COLLECTION AND HANDLING section and consult the specific TOSOH AIA System Operator's Manual, respectively. Barcoded primary tubes as well as sample cups can be run on AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800, AIA-2000 and AIA-360. When diluted samples are used on the AIA-360, the samples must be manually diluted in sample cups and the sample cups should be placed on the AIA-360. The primary tubes should not be used for diluted samples on the AIA-360.

B) Assay Procedure

1. Ensure a sufficient quantity of ST AIA-PACK BMG test cups for the number of samples to be run.
2. Load patient samples as instructed in the Operator's Manual and proceed with analysis. Note: The AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800 and AIA-2000 will require AIA-PACK SAMPLE TREATMENT CUPS if onboard dilutions are utilized.

PROCEDURAL NOTES

1. Lyophilized substrate must be completely dissolved.
2. Ligand assays performed by the TOSOH AIA System Analyzers require that the laboratory use water designated by the College of American Pathologists as Class I or by NCCLS as Type I. Water should be tested at least once per month and should be free of particulate matter including bacteria. The pH of the water should also be routinely tested. For further information, consult the NCCLS document "Preparation and Testing of Reagent Water in the Clinical Laboratory," NCCLS Document C3-A3, Volume 11 No. 13, originally approved as a guideline by NCCLS in October 1997.
3. If a specimen β_2 microglobulin concentration is found to be greater than 20 mg/L in serum or heparinized plasma specimens or 2.0 mg/L in urine specimens, the specimen should be diluted with the BMG Sample Diluting Solution and reassayed according to the Assay Procedure. The recommended dilution for specimens containing greater than 20 mg/L (serum or heparinized plasma) or 2.0 mg/L (urine) is 1:10. It is desirable to dilute the specimen so that the diluted specimen reads between 2.5 and 15 mg/L in serum and heparinized plasma specimens, and between 0.25 and 1.5 mg/L in urine specimens. If the autodilution feature is not being used, 51-fold or 5-fold dilution must be made for serum or heparinized plasma specimens and for urine specimens, respectively. For further information on the dilution of specimens, refer to the TOSOH AIA System Operator's Manual.
4. The TOSOH AIA System Analyzers can store two different calibration curves for each analyte at one time. Therefore, up to two different lots of ST AIA-PACK BMG test cups can be used during the same run.
5. If the assay specifications for this test are not ready in the system software, the specifications must be entered under test code 065.

CALCULATION OF RESULTS

The TOSOH AIA System Analyzers perform all sample and reagent handling operations automatically. The TOSOH AIA System Analyzers read the rate of fluorescence produced by the reaction and automatically convert the rate to β_2 microglobulin concentration in mg/L.

For samples requiring dilution, the AIA Nex-IA / AIA-21, AIA-600 II, AIA-1800 and AIA-2000 will automatically perform dilution and calculate results if the dilution factors are entered into the software. Dilution factors may be entered into the Test File, or pre-defined dilution factors may be selected in Specimen Processing.

EVALUATION OF RESULTS

Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially available control samples should be assayed according to the local regulations.

The minimum recommendations for the frequency of running internal control material are:

After calibration, three levels of the internal control are run in order to accept the calibration curve.

The three levels of controls are repeated when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustment or change).

After daily maintenance, at least two levels of the control should be run in order to verify the overall performance of the TOSOH AIA System Analyzers.

If one or more control sample value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the strict regulatory agency under which the laboratory operates.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, therapy, etc.).
- Using ST AIA-PACK BMG, the highest concentration of β_2 microglobulin measurable by 51-fold dilution (for serum or heparinized plasma specimens) or 5-fold dilution (for urine specimens) with the BMG Sample Diluting Solution is 0.4 mg/L, which corresponds to 20 mg/L in serum and heparinized plasma, and 2.0 mg/L in urine. The lowest measurable concentration of β_2 microglobulin is 0.002 mg/L (assay sensitivity), which corresponds to 0.10 mg/L in serum and heparinized plasma, and 0.010 mg/L in urine. Specimens that read < 0.002 mg/L may be repeated on a lesser dilution, and results calculated using the used dilution factor.
- Although the approximate value of the highest calibrator is 0.42 mg/L, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 0.4 mg/L, which corresponds to 20 mg/L in original serum and heparinized plasma specimens (dilution factor 51), and 2.0 mg/L in original urine specimens (dilution factor 5).
- Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- Lipemia has an insignificant effect on the assay except in the case of gross lipemia where spatial interference may occur.
- Specimens from patients taking medicines and/or medical treatment may show erroneous results.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show falsely elevated values when tested for β_2 microglobulin.
- For a more complete understanding of the limitations of this procedure, please refer to the SPECIMEN COLLECTION AND HANDLING, WARNINGS AND PRECAUTIONS, STORAGE AND STABILITY, and PROCEDURAL NOTES sections in this insert sheet.

EXPECTED VALUES

Each laboratory should determine a reference interval which corresponds to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted with regard to concomitant medications administered to the patient (26).

Reference Ranges

The intervals given here were determined in 161 serum samples and 54 urine samples from apparently healthy Asian individuals.

Reference Interval	
Serum	0.85 - 1.62 mg/L
Urine	0.027 - 0.265 mg/L

PERFORMANCE CHARACTERISTICS

ACCURACY

- a. Recovery: Three serum pools and three urine pools were spiked with three different levels of β_2 microglobulin and assayed before and after spiking.

Sample	Initial Value (mg/L)	BMG Added (mg/L)	Expected Value (mg/L)	Measured Value (mg/L)	Percent Recovery (%)
Serum A1	0.021	0.094	0.115	0.116	100.9
	0.021	0.047	0.068	0.065	95.6
	0.021	0.023	0.044	0.042	95.5
Serum B1	0.021	0.094	0.115	0.120	104.3
	0.021	0.047	0.068	0.068	100.0
	0.021	0.023	0.044	0.044	100.0
Serum C1	0.026	0.094	0.120	0.127	105.8
	0.026	0.047	0.073	0.072	98.6
	0.026	0.023	0.049	0.050	102.0
Urine A1	0.021	0.248	0.269	0.256	95.2
	0.021	0.124	0.145	0.138	95.2
	0.021	0.062	0.083	0.076	91.6
Urine B1	0.000	0.248	0.248	0.240	96.8
	0.000	0.124	0.124	0.110	88.7
	0.000	0.062	0.062	0.054	87.1
Urine C1	0.015	0.248	0.263	0.250	95.1
	0.015	0.124	0.139	0.129	92.8
	0.015	0.062	0.077	0.069	89.6

- b. Dilution: Three serum samples and three urine samples containing high concentrations of β_2 microglobulin were serially diluted with the BMG Sample Diluting Solution and assayed.

Sample	Dilution Factor	Expected Value (mg/L)	Measured Value (mg/L)	Percent Recovery (%)
Serum A2	none		0.397	
	7.5/10	0.298	0.285	95.6
	5.0/10	0.199	0.188	94.5
	2.5/10	0.099	0.091	91.9
	1.0/10	0.040	0.036	90.0
Serum B2	none		0.352	
	7.5/10	0.264	0.256	97.0
	5.0/10	0.176	0.162	92.0
	2.5/10	0.088	0.079	89.8
	1.0/10	0.035	0.031	88.6
Serum C2	none		0.340	
	7.5/10	0.255	0.261	102.4
	5.0/10	0.170	0.178	104.7
	2.5/10	0.085	0.087	102.4
	1.0/10	0.034	0.036	105.9
Urine A2	none		0.377	
	7.5/10	0.283	0.272	96.1
	5.0/10	0.189	0.185	97.9
	2.5/10	0.094	0.086	91.5
	1.0/10	0.038	0.035	92.1
Urine B2	none		0.383	
	7.5/10	0.287	0.269	93.7
	5.0/10	0.191	0.182	95.3
	2.5/10	0.096	0.086	89.6
	1.0/10	0.038	0.035	92.1

Urine C2	none		0.413	
	7.5/10	0.310	0.303	97.7
	5.0/10	0.207	0.205	99.0
	2.5/10	0.103	0.097	94.2
	1.0/10	0.041	0.038	92.7

PRECISION

- a. Within run precision was determined using three serum controls and two urine controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean (mg/L)	Pooled SD (mg/L)	CV (%)
Serum A3	0.023	0.00074	3.2
Serum B3	0.103	0.00250	2.4
Serum C3	0.247	0.00504	2.0
Urine A3	0.028	0.00047	1.7
Urine B3	0.227	0.00397	1.7

- b. Total precision was determined by the duplicate assay of three serum controls and two urine controls in 20 separate runs. The means of each run were used to calculate the pooled standard deviation (SD) and coefficient of variation (CV).

Sample	Mean (mg/L)	Pooled SD (mg/L)	CV (%)
Serum A3	0.023	0.00109	4.7
Serum B3	0.103	0.00386	3.7
Serum C3	0.247	0.00932	3.8
Urine A3	0.028	0.00087	3.1
Urine B3	0.227	0.00745	3.3

CORRELATION

The correlation between serum (x) and heparinized plasma (y) on ST AIA-PACK BMG was carried out using 67 patient specimens.

Slope	0.982
y-Intercept	-0.096
Correlation Coefficient	0.986
Range of Samples	0.71-19.97
Number of Samples	67

SENSITIVITY

The minimal detectable concentration (MDC) of β_2 microglobulin is estimated to be 0.002 mg/L which corresponds to 0.10 mg/L in original serum and heparinized plasma specimens (dilution factor 51) and corresponds to 0.01 mg/L in original urine specimens (dilution factor 5). The MDC is defined as the concentration of β_2 microglobulin which corresponds to the rate of fluorescence that is two standard deviations from the mean rate of fluorescence of 20 replicate determinations of a CALIBRATOR (1).

INTERFERENCE

Interference is defined, for the purposes of this study, as recovery outside of 10% of the known mean concentration of the specimen after the following substances are added to human specimens.

- Hemoglobin (up to 390 mg/dL), free bilirubin (up to 17 mg/dL) and conjugated bilirubin (up to 18 mg/dL) do not interfere with the assay.
- Lipemia, as indicated by triglyceride concentration (up to 1,600 mg/dL), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dL) does not interfere with the assay.
- Protein, as indicated by human albumin concentration (up to 5.0 g/dL), does not interfere with the assay.
- Heparin (up to 100 U/mL) does not interfere with the assay.

REFERENCES

- Berggard, I., and Bearn, A. G., Isolation and Properties of a Low Molecular Weight β_2 -Microglobulin Occurring in Human Biological Fluids. J. Biol. Chem. 243:4095 (1968).
- Nakamoro, K., et al., Multiple Common Properties of β_2 -Microglobulin and the Common Portion Fragment Derived from HLA Antigen Molecules. Proc Natl Acad Sci USA 70:2863-2865 (1973).
- Grey, H. M., Kubo, R. T., et al., The Small Subunit of HLA is β_2 -Microglobulin. J. Exp. Med. 138:1608 (1973).
- Peterson, P. A., et al., Proc Natl Acad Sci USA 71:35-39 (1974).
- Ervin, P. E., Wibell, L., Scand. J. Clin. Lab. Invest. 29:69-73 (1972).
- Schentag, J. J., et al., Early Detection of Aminoglycoside Nephrotoxicity with Urinary Beta-2-Microglobulin. J. Med. 9:201-210 (1978).
- Fleming, J. J., et al., Increased Urinary β_2 -Microglobulin after Cancer Chemotherapy. Cancer Treat. Rep. 64:581-588 (1980).
- Merle, L. J., et al., Renal Injury in Patients with Rheumatoid Arthritis Treated with Gold. Clin. Pharmacol. Ther. 28:216-222 (1980).
- Kjellstrom, T., Piscator, M., Phadedoc Diagnostic Communications No.1 (1979) Pharmacia Diagnostics AB, Uppsala, Sweden.
- Wibell, L., et al., Noninvasive Monitoring of Renal Transplant Function by Analysis of beta2-Microglobulin. Nephron 10:320-331 (1973).
- Shea, P. H., et al., Nephron 29:30-35 (1981).
- Bernier, G. M., Post, R. S., β_2 -Microglobulin. A Marker of Renal Homograft Survival, Transplantation 15:176 (1973).
- Edwards, L. C., et al., Noninvasive Monitoring of Renal Transplant Function by Analysis of beta 2-Microglobulin, Kidney International 23:767 (1983).
- Morrell, A., Riesen, W., Serum β_2 -Microglobulin, Serum Creatinine, and Bone Marrow Plasma Cells in Benign and Malignant Monoclonal Gammopathy, Acta Haematol (Basel) 65:87 (1980).
- Bataille, R., et al., Prognostic Factors and Staging in Multiple Myeloma: A Reappraisal. J. of Clin. Oncol. 4:80 (1986).
- Cassuto, J. P., et al., β_2 -microglobulin, a Tumour Marker of Lymphoproliferative Disorders. Lancet 2:950 (1978).
- Simonsson, B., β_2 -Microglobulin in Chronic Lymphocytic Leukemia, Scand. J. Haematol. 24:174 (1980).
- Crisp, A. J., et al., β_2 -Microglobulin Plasma Levels Reflect Disease Activity in Rheumatoid Arthritis. J. of Rheumatology 10:954 (1983).
- Weissel, M., et al., Serum β_2 -Microglobulin and SLE. Arthritis Rheum. 19:968 (1976).
- Strom, T., Serum β_2 -Microglobulin in Sjogren's Syndrome, Scand. J. Rheumatol. 7:97 (1978).
- Descos, L., et al., Serum levels of β_2 -Microglobulin - a New Marker of Activity in Crohn's Disease. New England J. Med. 301:440 (1979).
- Bhalla, R. R., Schwartz, M. K., Routine use of β_2 -Microglobulin Assays, Clin. Chem. 33:1957 (1987).
- Lambin, P., et al., Neopterin and β_2 -Microglobulin in Serum of HIV Seropositive Subjects During a Two-Year Follow-up. Clin. Chem. 34:1367 (1988).
- Fahey, J. L., et al., Enhanced Prediction of AIDS in HIV Seropositive Asymptomatic Men: Relative Value of 8 Cellular and Serologic Markers Singly and in Combination. 5th Int. AIDS Conference, Montreal, W. A. P. 62 (June 1989).
- Edelson, R., et al., Comparison of Activation Markers for Progression to AIDS, Ibid., W.A.P.76.
- Young, D., Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, Washington, DC, American Association for Clinical Chemistry Press (1990).

TOSOH CORPORATION
Shiba-koen First Bldg.
3-8-2 Shiba, Minato-ku, Tokyo 105-8623
Japan
Phone :+81-3-5427-5181
Fax :+81-3-5427-5220

TOSOH EUROPE N.V.
Transportstraat 4
B-3980 Tessenderlo, Belgium
Phone :+32-13-66 88 30
Fax :+32-13-66 47 49

TOSOH BIOSCIENCE, INC.
6000 Shoreline Ct., Suite 101
South San Francisco, CA94080, USA
Phone :+1-650-615-4970
Fax :+1-650-615-0415
Phone : (800)248-6764
Fax : (800)685-7595